Introduction

MStat 2 is has an entirely new patient evacuation module for nursing homes and other non-hospital facilities. The process is simple:


- 1. Get your patient or resident list into a simple Excel template that can be downloaded from MStat
- 2. Ensure that you have the basic required data elements and upload the spreadsheet into MStat.
- 3. You have options to update or replace patient/resident data already in MStat
- 4. You can manage groups of persons at once using a powerful Bulk Edit function

The Patient Template



Once logged in to EMSTAT, click the "Patient" link on the main page.


Elements of the Patient View



The basic elements of the Patient Data view include:

- 1. Along the top of the page, the view lists summary numbers of patients/residents by disposition
- The highlighted tab is the current view showing the basic data elements for each person- Patient
- 3. The Patients Expanded view shows additional data elements for each person
- 4. Click the Download Template to get an empty Excel template
- 5. Once you have your residents in the template, you can upload them in this view
- 6. The Medical Equipment view lets you indicate whether your resident is on a ventilator or other medical devices
- 7. The Evacuation tab lets you "move" residents to various locations during an evacuation
- 8. The Bulk Update function lets you perform critical actions such as changing locations on groups of residents.
- 9. Finally, a facility can add a patients one at a time if needed.

The Download Template



Clicking the "Download Template" button automatically downloads a blank template to your local computer. The template is named "NH_Patients_Upload-Template.xls"

The template data elements are listed below:

Field Name	Mandatory	Allowed Values		
LastName	Yes	text		
FirstName	Yes	text		
MiddleName	No	text		
MRNumber	Yes	text		
Gender	Yes	(Unknown,Male,Female)		
DOB	No	Date		
VentilatorDependent	No	(1,0) 1= Yes; 0 = No		
OxygenDependent	No	(1,0) 1= Yes; 0 = No		
OtherElectricalDeviceDependent	No	(1,0) 1= Yes; 0 = No		
Disposition	No	(In Facility, On Leave, Deceased, In Hospital, Evacuated)		
EvacueeType	No	(Staff,Patient,Guest,Shelteree)		
TransportationTraiage	No	(Red, Yellow, Green)		


Note there are 4 mandatory fields for a successful upload.

Entering Sample Data

	Α	В	С	D	E	F	G	
1	LastName	FirstName	MiddleName	MRNumber	Gender	DOB	VentilatorDependent	Охуд
2	text	text	text	text	('Unknown','Male','Female')	date	('1','0')	('1','
3	Yennie	Henry		444333aa	Male			
4	Prats	Rosanne		45454ffdd	Female			
5	Tietje	Malcolm		5675BAB	Male			
6								
7		Only mandatory data is included here						
8								
9								
10			×.					


In this example, we have populated a few persons in the template.

Uploading the Spreadsheet



In this step, we take the spreadsheet prepared in the last step and upload it to the system. We start by clicking the "Upload Data" tab.

The Upload Data Page




On this page, we start the process by choosing from one of three options:

- 1. Completely replace the list currently in the system
- 2. Add or update the list with patients in the spreadsheet
- 3. Delete the patients currently in the system.

In this example, we will choose "Add/Update Patients"

We will then click "Select" to find the template file with our patient data.

Uploading the Patient File



Once we locate the file, we click the "Upload Files" button.


Evaluating the Results



If you are successful, you will get a view similar to that shown above.


If there were errors, you would get a count of the "Total No of Records Failed" along with a listing of the errors.

Viewing our Patients



Once the upload is complete, you can view the imported patients by clicking on one of the view tabs. The view shown above is the "Patients" view


Bulk Update Function



In this example, we will select two patients to update. We first select them by clicking the check box next to their names and then clicking the "Bulk Update" button.


Bulk Update Actions

Bulk Update - Medical Evacuation Data




In this view, you can alter any one of the appropriate fields. Here we select Disposition and make a change.

Making the Change



Once a selection is made, click "Update".

Viewing The Bulk Update Results



If we go back to "Patient", we can see the result of our Bulk Update.

Viewing All Persons



After clicking the "View All" link, we can see all persons in our system. We can now perform other bulk update functions as dictated by our evacuation progress.

Also note that at any time you can export your patients and their data to Excel.

Conclusion

This has been a high level overview of the Patient management function. As you experiment with this feature, you will notice error messages if you don't make proper choices in the process.

We will add details to this guide over the coming months, so stay tuned!